



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	904C
GENERAL DESCRIPTION	SRK 10ML 16X100 LARGE PL.250PK
DETAILED DESCRIPTION	10 boxes each containing 25 paper-film pouches each containing one tube with green screw cap and 10 ml of solution SRK and at the side a plastic applicator with rayon tip.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
SRK	Sodium Chloride Potassium Chloride Calcium Chloride Lecithin Sodium Thiosulphate Sodium Thioglycollate Sodium Pyruvate Sodium Bicarbonate Polysorbate 80

SAMPLING TUBE	TUBE WITH ROUNDED BOTTOM	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 101.80 mm \pm 0.20 mm TUBE DIAMETER: 16.10 mm \pm 0.20 mm INTERNAL DIAMETER TUBE: 14.10 mm \pm 0.20 mm MATERIAL: POLYPROPYLENE
CAP	GREEN CAP	DIMENSIONAL SPECIFICATIONS
		LOWER CAP DIAMETER: 21.50 mm \pm 0.20 mm CAP LENGTH: 16.30 mm \pm 0.20 mm HIGHER CAP DIAMETER: 20.50 mm \pm 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.02-30/04/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	PLASTIC APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STEM LENGTH: 190.0 mm \pm 1.0 mm BREAK POINT: 116.50 mm \pm 2.50 mm DIAMETER HANDLE PART: 4.00 mm \pm 0.08 mm TOTAL LENGTH: 194.50 mm \pm 5.00 mm LENGHT OF THE FIBRE TIP: 24.0 mm \pm 3.0 mm FIBRE TIP DIAMETER: 10.00 mm \pm 0.50 mm

PACKAGING SPECIFICATION

INNER BOX QTY	25
OUTER BOX QTY	250
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	24 MONTHS
MEDIUM pH	7.0 \pm 0.5