

Lab-Aire® II Drying Racks

Laboratory space is always in high demand making the small footprint of the Lab-Aire® II Drying Racks a practical solution for drying and storing fragile, frequently used labware.

Choose from **Non-electric Benchtop** or **Wallmount** drying racks or accelerate drying with the warm, flowing air from **Electric Benchtop** models. All styles have removable pegs in 3 different lengths that can be positioned for optimal use. Extra accessories and a variety of different sizes/configurations let you create a drying system that works for you.

Wallmount - Non-electric

- Space-saving, move drying off the bench
- Easily mounts on walls or standard pegboard
- Single-sided; can be custom made with additional height

Benchtop - Non-electric

- Portable, use anywhere convenience
- Offers the most choices in size
- Single- and Double-sided models

Benchtop - Electric

- Accelerate drying with warm, circulating air
- Use with or without power cycle
- Single- and Double-sided models: 120V or 230V



Wallmount - Non-electric
Single-sided, 2 tier



Benchtop - Non-electric
Double-sided, 1 tier



Benchtop - Electric
Single-sided, 3 tier

Increase Drying Space Three Ways

1. Go up in height – 1, 2 and 3 tier models available in all styles
2. Double-sided benchtop models double your capacity with minimal increase in footprint
3. Use multiple Lab Aire® II Dryers

Increase Flexibility

Each Lab Aire® II Dryer comes with an assortment of pegs in various lengths. Electric models also have baskets for holding small labware. Additional pegs and baskets* to customize your dryer can be purchased separately.

* Baskets fit electric models only.

Lab-Aire® II Dryers Quick Comparison Chart

Select from 12 different models with up to 114 peg positions.

- Rugged ABS plastic will never corrode and easily wipes clean
- Drip/storage trays with scalloped wells hold long thin labware
- Small (3"), Medium (4"), and Large (5") pegs can be moved and repositioned to hold a variety of labware including beakers, flasks, bottles, test tubes and cylinders
- Extra pegs and drying baskets available

| | Style | Wallmount | | | Dimensions (W x D x H) | Peg Slots | Peg Sizes* | Model |
|---------------------------|-------------------|-----------|----------|----------|--|--------------|--|--|
| | | | Benchtop | Electric | | | | |
| Wallmount Non-Electric | Single 1 Tier | X | X | | 37.5 x 12.7 x 24cm (14.75 x 5 x 9.4") | 19 | S - 6, M - 7, L - 6 | 18933-0011 |
| | Single 2 Tiers | X | | | 37.5 x 12.7 x 41.7cm (14.75 x 5 x 16.4") | 38 | S - 12, M - 14, L - 12 | 18933-0014 |
| | Single 3 Tiers | X | | | 37.5 x 12.7 x 59.4cm (14.75 x 5 x 23.4") | 57 | S - 18, M - 21, L - 18 | 18933-0013 |
| Benchtop Non-Electric | Single 1 Tier | X | X | | 37.5 x 12.7 x 24cm (14.75 x 5 x 9.4") | 19 | S - 6, M - 7, L - 6 | 18933-0011 |
| | Single 2 Tiers | | X | | 37.5 x 17.8 x 39cm (14.75 x 7 x 15.4") | 38 | S - 12, M - 14, L - 12 | 18933-0012 |
| | Single 3 Tiers | | X | | 37.5 x 17.8 x 57cm (14.75 x 7 x 22.4") | 57 | S - 18, M - 21, L - 18 | 18933-0015 |
| | Double 1 Tier | | X | | 37.5 x 25.4 x 21.3cm (14.75 x 10 x 8.4") | 38 | S - 12, M - 14, L - 12 | 18933-0021 |
| | Double 2 Tiers | | X | | 37.5 x 25.4 x 39cm (14.75 x 10 x 15.4") | 76 | S - 24, M - 28, L - 24 | 18933-0022 |
| | Double 3 Tiers | | X | | 37.5 x 25.4 x 57cm (14.75 x 10 x 22.4") | 114 | S - 36, M - 42, L - 36 | 18933-0023 |
| Benchtop Electric | Single 2 Tiers | X | | X | 42.5 x 19 x 40cm (16.75 x 7.5 x 15.7") | 19 | S - 6, M - 7 L - 6 2 Baskets, 7 Peg Plugs | 18819-0012 (120V) 18819-2301 (230V) |
| | Single 3 Tiers | X | | X | 42.5 x 19 x 57.7cm (16.75 x 7.5 x 22.7") | 38 | S - 12, M - 14 L - 12 2 Baskets, 14 Peg Plugs | 18819-0013 (120V) 18819-2302 (230V) |
| | Double 2 Tiers | X | | X | 42.5 x 25.4 x 40cm (16.75 x 10 x 15.7") | 38 | S - 12, M - 14 L - 12 4 Baskets, 14 Peg Plugs | 18819-0022 (120V) 18819-2303 (230V) |
| | Double 3 Tiers | X | | X | 42.5 x 25.4 x 57.7cm (16.75 x 10 x 22.7") | 76 | S - 24, M - 28 L - 24 4 Baskets, 28 Peg Plugs | 18819-0023 (120V) 18819-2304 (230V) |

*S - Small • M - Medium • L - Large