



MINI MICROCENTRIFUGE MULTISPIN 12



MULTISPIN 12		
External dimensions		
Width	mm	176
Height	mm	122
Depth	mm	196
Weight	kg	1,5
Performance		
Max speed	rpm	12.000
Speed increments	rpm	1.000
Max RCF	xg	7.900
Timer	min.	1÷99
Max capacity	ml	6 x 2 ; 6 x 0,5 ; 6 x 0,2
Rotor heads	type	A 6-2 (6 x 2 ; 6 x 0,5 ; 6 x 0,2 ml)
	type	PCR 2-8 (2 strips of 8 PCR tubes of 0,2 ml)
Voltage	Volt	230
Noise level	db	≤ 45db