



## Quick Overview

Removes organic - particulates and chlorine. contaminants to protect RO module. 2 cartridges are filters 5 $\mu$  and 2 cartridges are activated carbon

<b>Part Number</b>	PCK200
<b>Compatible with</b>	Micromatic (2 changes)
<b>Part Number</b>	PCK200
<b>Unit of measure</b>	Individual

Distribué par :

Z.A de Gesvrine - 4 rue Képler - B.P.4125  
44241 La Chapelle-sur-Erdre Cedex - France  
t. : +33 (0)2 40 93 53 53 | f. : +33 (0)2 40 93 41 00  
commercial@humeau.com



w w w . h u m e a u . c o m