



# Product Specifications

**CL01.1105**

IONEX

## **Potassium standard solution 1000 µg/ml for IC (Spectro ECON)**

Potasio solución patrón  
Potassium solution étalon  
Kalium standaard oplossing  
Kalium Standard Lösung  
Potassio soluzione standard

**Solution contains 990 - 1010 µg/ml of K - Method Accuracy ±0.5%**  
**1.907 g KCl / 1 H<sub>2</sub>O**  
**Traceability: Reference standards acc. NIST SRM 3141a**

Mol weight: 74.56 g/mol  
Density: 1,00 g/ml  
HS Nr: 38220000  
CAS Nr: 7447-40-7  
UN:  
ADR:

### **Specifications:**