

# 2- glass door tall safety cabinet with filtration - 2 compartments

## COMPLIANCE

NF X 15-211, Electric fan in compliance with CE specifications

## ADVANTAGES

- CE certified
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics

### BUT ALSO...

- Lighting control allowing to survey the proper working of the fan
- Front control window to see immediately the filter
- Stop 99,97% of the particles with a diameter bigger than 0.3 µm
- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Normalized labelling provided with the equipment
- Available in 2 models : glass doors or full doors
- Fan in compliance with CE specifications
- Active coal filters adsorption capacity and filtration system certified NF X 15-211
- No modification of building's structure required
- Easy and aesthetic installation
- Active coal filter tested to the standard NF X 15-211
- Easy replacement of the active coal filter



## TECHNICAL FEATURES

Model	<b>AF5-2X</b>	Door closure	<b>manual closing</b>
Compliance	<b>NF X 15-211, Electric fan in compliance with CE specifications</b>	Shipping dimensions HxWxD (mm)	<b>2340 x 1200 x 800</b>
Cabinet materials	<b>10/10th epoxy painted steel structure</b>	Ventilation	<b>Inlet located at the bottom of the cabinet, outlet located at the top</b>
Doors color	<b>white RAL 9010</b>	Number of compartments	<b>2</b>
Number of doors	<b>2</b>	Containment sump	<b>2</b>
Doors type	<b>swing doors</b>	HS Code	<b>94032080</b>
Hinges	<b>right &amp; left</b>		

## COMPLEMENTARY PRODUCTS

- Activated carbon filter for ammonia vapors - for filtration box H50C (or old cabinets range 12 before 2002) (Réf : AMM50)
- Activated carbon filter for organic and corrosive vapors - for filtration box H50C (or old cabinets range 12.X before 2002) (Réf : CORG51)
- Activated carbon filter for formaldehyde vapors - for filtration box H50C (or old cabinets range 12.X before 2002) (Réf : FOR50)
- Activated carbon filter for organic vapors - for filtration box H50C (or old cabinets range 12.X before 2002) (Réf : ORG50)

## OPTIONS

- Retention shelf (Réf : EE1)
- Containment sump in PVC (Réf : SBPVC)
- Shelf protection in PVC (Réf : SEPVC)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Weight (kg)
AF5-2X	2135 x 1110 x 545	(2 x )1795 x 555 x 480	130.00