

XS BALANCE mod. BL 6001

Basic -6000gr. - 0,1 Code: 27000293

Réf.: 042.270002.93



DATA SHEET

Regulations

Regulations CE

Performance

| | |
|----------------------------|---------|
| Capacity | 6000 g |
| Divisions | 0,1 g |
| Linearity | ± 0,2 g |
| Min. weight | 1 g |
| Response time | 1,5 sec |
| Ripetibility | 0,2 g |
| Working temperature (max.) | 35 °C |
| Working temperature (min.) | 5 °C |

Feature

| | |
|---------------------------|---------------|
| Calibration | External |
| Display | Backlight LCD |
| Ripetibility | 0,2 g |
| Selectable measuring unit | Yes |
| Sensor | Nuosheng |

Supply

| | |
|------------------------|---------------|
| Power supply | 220 V - 50 Hz |
| Rechargeable batteries | Yes |

Dimension and Weight

| | |
|----------------------------|--------------|
| External dimensions | 280x180 mm |
| Net weight | 1,2 Kg |
| Working plate - dimensions | 143 x 192 mm |