


Specification

Incuwater-Clean™

A5219

Product Code:	A5219
Product Name:	Incuwater-Clean™
Headline Comment:	<ul style="list-style-type: none"> • Incuwater-Clean™ is a registered trademark of AppliChem GmbH.
Specifications:	<ul style="list-style-type: none"> : • for the prevention of microbial growth in incubator water baths. : • Disinfectant solution for CO₂-incubator water baths. : • 100X concentrated solution. Use 50 ml per 5 liter of incubator water bath.
WGK:	1
Storage:	2 - 8°C
CS:	38221900
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340