



Midisart® 2000

Ready-to-use filter for sterilizing gases and venting

Product Information

The Midisart® 2000 is the best choice for small-scale sterile air | gas filtration applications. Superior flow rates and seamless integration into single-use assemblies provide Midisart® customers with a distinct advantage that directly translates into cost savings. Individual 100% factory integrity testing verifies the quality of every filter prior to use. Moreover our new production-line with its enhanced capacity ensures a quick and reliable supply at any time.



Applications

Midisart® are excellent for

- Sterile venting of filling vessels and fermentation carboys, including culture vessels and CO₂ incubators (6 to 120 liters)
- Venting of holding tanks for sterile, distilled water and liquid culture media
- Autoclave venting
- In-line sterilization of and particulate removal from air and gases, such as sterilization of air for small fermenters

Flexibility

Midisart® 2000 were designed for maximum handling ease and safety. Tapered hose barbs ensure a simple and secure hold for 6- to 12-mm inner diameter tubing. Other connector types such as small hose barb (for tubings with 4–12 mm inner diameter), 8" NPT thread and TriClamp are also available. Midisart® is lightweight – only 20 g – so it will not weigh down or kink the tubing.

Reliability and safety

Midisart® 2000 is integrity testable and delivers reproducible results. The membrane is reinforced with polypropylene gauze, giving the Midisart® unit added stability and making it pressure resistant up to 3 bar | 44 psi. Midisart® 2000 entirely eliminates moisture breakthrough because of its inherently hydrophobic PTFE material. In addition, Midisart® is bio safe because all materi-

als of construction meet the requirements of the current USP Plastics Class VI testing. Midisart® 2000 units easily withstand at least 20 autoclaving cycles with no loss in performance.

Quality Control

Each Midisart® 2000 is tested by integrity prior to release. This 100% test has been defined as the final release criteria. The lot number and the individual unit number are lasered on the top part of each Midisart® 2000 housing to ensure complete traceability.

Performance

With a diameter measuring just 64 mm, Midisart® incorporates a filter area of 20 cm², which means that it is "packed" with high flow rate performance power. Midisart® multiply filtration performance in more ways than one.

Lot Number and Individual Unit Number

Lasered on the top part of each Midisart® - housing.

Integrity Testing

Midisart® 2000 filters can be integrity tested (Bubble Point Test, 60% IPA, BPmin. = 1.1 bar). The Midisart® Test Manifold (Order code 1Z-LB-0002) can be used to test 10 filters in parallel which simply leads to a reduction of test time by 90%.

Technical Data

Specifications

Materials

Membrane PTFE – reinforced with polypropylene gauze
Housing Polypropylene

Pore Sizes

0.2 μm

0.45 μm

Filtration Area

20 cm^2 | 3 square inch

Housing diameter

64 mm | 2,5"

Connectors

Hose Barb, 1/8" NPT, Triclamp, small Hose Barb

Biosafety

Materials pass current UPS Plastics, class VI
Non fiber releasing according to 21 CFR

Sterilization

Max. Temperature 134°C

Max. Autoclaving 60 Cycles

Not suitable for gamma irradiation

Gamma stable: see Midisart® BV

Technical References

Validation Guide

SLD5702-e

Ordering Information

Order No.	Pore Size [μm]	Membrane	Connectors	Qty.	Sterile
17804-----E	0.45	PTFE	Hose Barb Hose Barb	12	Yes
17804-----G	0.45	PTFE	Hose Barb Hose Barb	25	Yes
17804-----NPE	0.45	PTFE	1/8" 1/8" NPT	12	Yes
17804-----NPG	0.45	PTFE	1/8" 1/8" NPT	25	Yes
17805-----E	0.2	PTFE	Hose Barb Hose Barb	12	Yes
17805-----G	0.2	PTFE	Hose Barb Hose Barb	25	Yes
17805-----NPE	0.2	PTFE	1/8" 1/8" NPT	12	Yes
17805-----NPG	0.2	PTFE	1/8" 1/8" NPT	25	Yes
17805-----UPN	0.2	PTFE	Hose Barb Hose Barb	100	No
17805-----UPQ	0.2	PTFE	Hose Barb Hose Barb	500	No
17805-----UQN	0.2	PTFE	Hose Barb Hose Barb	100*	No
17809-----UNN	0.2	PTFE	1/8" 1/8" NPT	100	No
17812-----UNN	0.2	PTFE	1/8" Hose Barb	100	No
17805-----TCN	0.2	PTFE	TriClamp TriClamp	100	No
17877-----UPN	0.2	PTFE	Small Hose Barb Small Hose Barb	100	No

* = plastic bucket – the bucket with its smooth surface can be easily transferred into clean-room areas without risking particle release from the packaging material.



Standard Hose Barb



Small Hose Barb



1/8" NPT Thread



TriClamp

Sales and Service Contacts

For further contacts, visit www.sartorius-stedim.com

Europe

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

Sartorius Stedim Systems GmbH
Robert-Bosch-Strasse 5 – 7
34302 Guxhagen

Phone +49.5665.407.0
Fax +49.5665.407.2200

France

Sartorius Stedim FMT S.A.S.
ZI des Paluds
Avenue de Jouques – CS 91051
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI des Paluds
Avenue de Jouques – CS 71058
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.846545

Austria

Sartorius Stedim Austria GmbH
Modcenterstrasse 22
1030 Vienna

Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Stedim Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.80
Fax +32.2.756.06.81

Hungary

Sartorius Stedim Hungária Kft.
Kagyló u. 5
2092 Budakeszi

Phone +36.23.457.227
Fax +36.23.457.147

Italy

Sartorius Stedim Italy S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Stedim Netherlands B.V.

Phone +31.30.60.25.080
Fax +31.30.60.25.099

filtratie.nederland@sartorius-stedim.com

Poland

Sartorius Stedim Poland Sp. z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn

Phone +48.61.647.38.40
Fax +48.61.879.25.04

Russian Federation

LLC "Sartorius ICR" and LLC "Biohit"
Uralskaya str. 4, Lit. B
199155, Saint-Petersburg

Phone +7.812.327.5.327
Fax +7.812.327.5.323

Spain

Sartorius Stedim Spain, S.A.U.
Avda. de la Industria, 32
Edificio PAYMA
28108 Alcobendas (Madrid)

Phone +34.902.110.935
Fax +34.91.358.96.23

Switzerland

Sartorius Stedim Switzerland AG
Ringstrasse 24 a
8317 Tagelswangen

Phone +41.52.354.36.36
Fax +41.52.354.36.46

U.K.

Sartorius Stedim UK Ltd.
Longmead Business Centre
Blenheim Road, Epsom
Surrey KT19 9 QQ

Phone +44.1372.737159
Fax +44.1372.726171

Ukraine

LLC "Biohit"
Post Box 440 "B"
01001 Kiev, Ukraine

Phone +380.44.411.4918
Fax +380.50.623.3162

America

USA

Sartorius Stedim North America Inc.
5 Orville Drive, Suite 200
Bohemia, NY 11716

Toll-Free +1.800.368.7178
Fax +1.631.254.4253

Argentina

Sartorius Argentina S.A.
Int. A. Ávalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Avenida Senador Vergueiro 2962
São Bernardo do Campo
CEP 09600-000 - SP- Brasil

Phone +55.11.4362.8900
Fax + 55.11.4362.8901

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente
No. 149
Ciudad Satélite
53100, Estado de México
México

Phone +52.5555.62.1102
Fax +52.5555.62.2942

Asia | Pacific

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

China

Sartorius Stedim Biotech (Beijing) Co. Ltd.
No. 33 Yu'an Road
Airport Industrial Park Zone B
Shunyi District, Beijing 101300

Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim (Shanghai)
Trading Co. Ltd
3rd Floor, North Wing, Tower 1
No. 4560 Jinke Road
Zhangjiang Hi-Tech Park
Pudong District
Shanghai 201210, China

Phone +86.21.68782300
Fax +86.21.68782332 | 68782882

Sartorius Stedim Biotech (Beijing) Co. Ltd.
Guangzhou Representative Office
Unit K, Building 23
Huihua Commerce & Trade Building
No. 80 Xianlie Middle Road
Guangzhou 510070

Phone +86.20.37618687 | 37618651
Fax +86.20.37619051

India

Sartorius Stedim India Pvt. Ltd.
#69/2-69/3, NH 48, Jakkasandra
Nelamangala Tq
562 123 Bangalore, India

Phone +91.80.4350.5250
Fax +91.80.4350.5253

Japan

Sartorius Stedim Japan K.K.
4th Fl., Daiwa Shinagawa North Bldg.
8-11, Kita-Shinagawa 1-chome
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.4331.4300
Fax +81.3.4331.4301

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park II
Singapore 117528

Phone +65.6872.3966
Fax +65.6778.2494

South Korea

Sartorius Korea Biotech Co., Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, BunDang-Gu
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700
Fax +82.31.622.5799



▶ www.sartorius-stedim.com