


# LE427

For semi-solid and solid samples as  
meat, fruits, cheese and soil

pH Ref °C ISM



<b>Sensor type</b>	combined pH electrode
<b>Order number</b>	51340333
<b>Specifications</b>	
pH range	1...11
zero point (pH)	7.00 ± 0.25
temperature range	0...50 °C
temperature probe	-
shaft material	glass
membrane glass	LoT
membrane resistance	< 250MΩ
reference system	Ag / AgCl
reference diaphragm	ceramic
reference electrolyte	gel
bridge electrolyte	-
electrode head	-
cable length	1.0 m
connector	BNC
cell constant (cm <sup>-1</sup> )	-
cell type	-
indicating metal	-
<b>Storage</b>	
3M KCl (51350072)	