


LE621

Robust DO sensor for lab use

DO °C ISM



Sensor type	galvanic DO
Order number	51340403
Specifications	
measuring range	0...20 mg/L
zero point (pH)	-
temperature range	0...50 °C
temperature probe	NTC 30 kΩ
shaft material	PVC
membrane glass	-
membrane resistance	-
reference system	-
reference diaphragm	-
reference electrolyte	-
bridge electrolyte	-
electrode head	-
cable length	1.0 m
connector	BNC / Cinch
cell constant (cm ⁻¹)	-
cell type	-
indicating metal	-
Storage	
10% NaCl	