


# LE703 IP67

Robust conductivity probe for outdoor applications

Cond °C ISM



<b>Sensor type</b>	conductivity cell
<b>Order number</b>	30266974
<b>Specifications</b>	
measuring range	0.01...200 mS/cm
zero point (pH)	-
temperature range	0...80 °C
temperature probe	NTC 30 kΩ
shaft material	epoxy
membrane glass	-
membrane resistance	-
reference system	-
reference diaphragm	-
reference electrolyte	-
bridge electrolyte	-
electrode head	-
cable length	2.0 m
connector	LTW
cell constant (cm <sup>-1</sup> )	0.55
cell type	4 graphite poles
indicating metal	-
<b>Storage</b>	
Dry	