



# Product specification

**Product number** 7.85420  
**Packaging number** 7.85420.782  
**Product name** NBB-PCR enrichment medium

**Version** 3.0  
**Date of version** 11.12.2018  
**Page** 1/2

---

## Information about composition and property

### Appearance and consistency

clear to slightly cloudy, brown-orange, broth

### Odour and taste

unfit for consumption, therefore not applicable

### Descriptions

#### Ingredients

water, beer, sugar, peptone, yeast extract, growth factors, buffer, inhibitors

#### Application

For the detection of beer spoiling bacteria via PCR

#### Indication

Reseal opened ware tightly and use quickly.

Please refer to user manual.

Unfit for consumption.

pH-value: 5,7 (5,5 - 5,9)

### Other information

This product contains allergenes subject to labelling according to Regulation (EU) No 1169/2011 and amendments. For details see description.

NBB®-PCR Broth is free of DNA of beer spoiling bacteria as well as PCR inhibitors in order to avoid false positive or false negative results.

---

## Storage and packaging information

Packaging and content: Card box, (9 x 250 ml)

### Shelf life in the unopened original packaging

At 4-8 °C : 540 day(s)

### Storage and handling

Chilled storage.

Store in a dark place.

Do not freeze.

### Microbiology

Sterilized and ready for use.



## Product specification

**Product number** 7.85420  
**Packaging number** 7.85420.782  
**Product name** NBB-PCR enrichment medium

**Version** 3.0  
**Date of version** 11.12.2018  
**Page** 2/2

---

### Dangerous good and hazardous material information

No dangerous good.  
No hazardous material.

---

### Disclaimer

This product specification describes the composition and respectively the properties of our product to the best of our knowledge. It replaces all previous versions and is valid without signature.