TITRALAB AS1000 SERIES SAMPLE CHANGER



Applications

- Food and Beverage
- Drinking Water
- Industrial Water
- Wastewater

The AS1000 Sample Changer is the perfect laboratory assistant for automating analyses on sample series.

Titration automation of sample series can be easy. Titralab[®] AS1000 series sample changer is used in connection with AT1000 Series titrators and TM1000 Windows Software to automatically perform sample series analyses (TM1000 Software and connection cables are included in the standard delivery).

Run sample analyses in total autonomy, without user intervention

The rotating module turns the tray to change the sample position while the lifting module moves sensors and tubes up and down into each sample. Single-radius sample trays for up to 30 samples.

Complete automation to free up your time

AS1000 Sample Changer multi-parameter automation solutions release operators' time from cumbersome repetitive analyses for more valuable tasks that cannot be automated.



System setup illustration with AT1000 Series Titrator and AS1000 Sample Changer



-5 to 40 °C (23 to 104 °F) Safety IEC/EN 61010-1;

EMC IEC/EN 61326-1

Warranty 1 year

*Subject to change without notice.

Technical Data*

Dimensions (H x W x D) 300 mm x 395 mm x 535 mm Weight Approximately 22.05 lbs. (10 kg) Instrument input: 24 VDC, 60 VA;

Power Requirements (Voltage)

External plug-in power supply:

24 VDC, 1.33 A

50/60 Hz

Power Requirements

(Hz)

Interface RS232 serial Stirrer Magnetic

Material Enclosures Lifting module: Enamelled steel

and Polypropylene;

Rotating module: Polypropylene

Order Information

Sample Changer Models

AS1000.97.12150 Titralab® AS1000 Series Sample Changer, 12 positions, 50 mL & 150 mL beakers, for two probes

AS1000.97.20090 Titralab® AS1000 Series Sample Changer, 20 positions, 90 mL beakers, for one probe AS1000.97.30050 Titralab® AS1000 Series Sample Changer, 30 positions, 50 mL beakers, for one probe

The scope of supply of each AS1000 Sample Changer includes TM1000 PC Software and connection cables.

Accessories and Spare Parts

LZE108 Set of beakers, 50 mL, PP (10 pcs) LZE109 Set of beakers, 150 mL, PP (10 pcs) **LZE193** Set of beakers, 50 mL, PP (30 pcs) **LZE194** Set of beakers, 90 mL PP (20 pcs) **LZE217** Magnetic bars, 12 x 4.5 mm (30 pcs) **LZE136** Magnetic bars, 20 x 6 mm (10 pcs) **LZE201** Carousel, 20 positions, 90 mL beakers **LZE202** Carousel, 30 positions, 50 mL beakers

LZE203 Carousel, 12 positions, 50 mL & 150 mL beakers

LZE213 Beaker ring adapter, 50 mL (12 pcs)

HACH World Headquarters: Loveland, Colorado USA

United States: 800-227-4224 tel 970-669-2932 fax orders@hach.com Outside United States: 970-669-3050 tel 970-461-3939 fax int@hach.com

hach.com

Printed in U.S.A.

©Hach Company, 2015. All rights reserved. In the interest of improving and updating its equipment, Hach Company reserves the right to alter specifications to equipment at any time

