



designed for scientists



BMT-20-S-M-gamma Tube for grinding with stainless steel balls

/// Data Sheet

Components: Tube, stainless steel balls, cover

Volume: 2 - 15 ml

Application: Dry milling of dry and brittle samples (e.g. kaolin, gypsum, coloured pigments, tablets), cell maceration, processing of materials mixed with fluids

Further suitable for:

- additional media quantities may be added during operation
- sampling during operation

www.ika.com

Subject to technical changes





designed for scientists

- sample preparation under sterile conditions
- germ free conservation of samples (tissue, blood, etc.)

Packing unit: 20 pieces





designed for scientists

Technical Data

Volume range (H2O) [l]	.002 - .015
Tube diameter [mm]	29
Volume Tube [ml]	20
Dry Operation	yes
FDA conformity	yes
Material: Balls	stainless steel
Material: Membrane	PP
Material: Testing Vessel and Cap	PP
Material: Gyrate Stick	PEEK
Sterile	yes
Weight [kg]	0.484

