

**CERTIFICATE OF ANALYSIS**

**PRODUCT** CM1082B  
**VIOLET RED BILE GLUCOSE AGAR (ISO)**

**LOT NUMBER** 3505173

**EXPIRY DATE** 2027.06.30

**DATE OF MANUFACTURE** 2022.06.13

**Delivery/Customer information**

Date Printed  
 2022.07.05

Delivery No.

Customer  
 Customer Order Number

Physical Characteristics	Results	Specification
Appearance	Straw/pink powder	Straw/pink powder
Colour on reconstitution	Purple	Purple
pH (25°C)	7.4	7.2 - 7.6
Clarity	Clear	Clear

Microbiological Performance	Control cfu	Test cfu	Recovery Test %	Description
<b>Aerobic incubation at 37 ± 2°C for 24 ± 2 hours</b>				
<b>Pour plate technique</b>				
<i>Klebsiella pneumoniae</i> ATCC®29665	69	76	110	Purple/pink cols & halo
<i>Proteus mirabilis</i> ATCC®12453	96	102	106	Purple cols & halo
<i>Shigella sonnei</i> ATCC®25931	105	86	82	Irregular, purple/pink cols & halo
<i>Enterobacter aerogenes</i> ATCC®13048	66	62	94	Purple/pink cols & halo
<b>Surface plate technique</b>				
<i>Pseudomonas aeruginosa</i> ATCC®9027	93	86	92	Straw cols, no halo

**Testing has been performed in accordance with ISO11133:2014****Aerobic incubation at 37 ± 2°C for 24 ± 2 hours****Pour plate technique**

<i>E. coli</i> ATCC®8739 WDCM00012	56	48	86	Purple/pink cols & halo
<i>E. coli</i> ATCC®25922 WDCM00013	97	90	93	Purple/pink cols & halo
<i>S. typhimurium</i> ATCC®14028 WDCM00031	66	57	86	Purple/pink cols & halo
<i>S. enteritidis</i> ATCC®13076 WDCM00030	73	67	92	Purple/pink cols & halo

**Surface plate technique**

<i>E. faecalis</i> ATCC®19433 WDCM00009	>1E+04	0	-	No growth
<i>E. faecalis</i> ATCC®29212 WDCM00087	>1E+04	0	-	No growth

Control Medium: Tryptone Soya Agar

A satisfactory result is represented by recovery of positive strains equal to or greater than 50% of the control medium.

Refer to product specification for full details.



Tested by the Quality Control Laboratory  
 OXOID LIMITED

Wade Road, Basingstoke, Hampshire RG24 8PW. England  
 oxid@thermofisher.com www.oxid.com

FDA Reg No. 8010096



## CERTIFICATE OF ANALYSIS

**PRODUCT** CM1082B  
**VIOLET RED BILE GLUCOSE AGAR (ISO)**  
**LOT NUMBER** 3505173

---

The information given is believed to be correct, all results reported in this certificate relate only to the product and lot stated in this certificate of analysis. However, both the information and the product are offered without warranty for any application other than that specified in the current Oxoid Manual. This certificate shall not be reproduced except in full, without written approval of the Quality Control Laboratory, Oxoid Limited, Basingstoke.

The results reported were obtained at the time of release.

Lot Accepted. 2022.07.04

.....  
Carissa Courtney  
Director Quality, EMEA

ATCC is a registered trade mark of the American Type Culture Collection.

NCTC and National Collection of Type Cultures are registered trade marks of the Health Protection Agency.



**Tested by the Quality Control Laboratory**  
**OXOID LIMITED**  
**Wade Road, Basingstoke, Hampshire RG24 8PW. England**  
**oxid@thermofisher.com www.oxid.com**